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ab Seite 3

## **Instruction manual**

starting on page 8

## **Mode d'emploi**

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## **Handleiding**

vanaf pagina 16

DE

GB

FR

NL

08916\_V0

## Inhalt

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## Contents

Meaning of the Symbols in These Instructions	9
Intended Use	9
Safety Instructions	9
Dangers of Injury	9
Notes For Your Health	9
Instructions for Use	10
Items Supplied	10
Use	10
Disposal	11

Dear Customer,

We are delighted that you have chosen our teeth whitening strips.

Before using the product for the first time, please read the instructions through carefully and keep them in a safe place. If the product is given to someone else, it should always be accompanied by these instructions. The manufacturer and importer do not accept any liability if the information in these instructions is not complied with!

We reserve the right to modify the product, packaging or enclosed documentation at any time in conjunction with our policy of continuous development.

If you have any questions, contact the customer service department via our website:

**[www.service-shopping.de](http://www.service-shopping.de)**

## Meaning of the Symbols in These Instructions

All safety instructions are identified with this symbol. Please read them carefully and comply with the safety instructions in order to prevent personal injury and damage to property.

Tips and recommendations are identified with this symbol.

## Intended Use

- ☐ This product is suitable for use on healthy teeth.
- ☐ The product is **not** suitable for children and animals. Keep out of the reach of children and animals.
- ☐ The product is **not** suitable for women who are pregnant or breastfeeding.
- ☐ Not suitable for allergy sufferers.
- ☐ The product must not be used on false teeth.
- ☐ The product should only be used as described in the instructions on the packaging. Any other use is deemed to be improper.
- ☐ The product is for personal use only and is not intended for commercial applications.
- ☐ The warranty does not cover any defects caused by incorrect use.

## Safety Instructions

### Dangers of Injury

- ☐ Keep children and animals away from the product and packaging!
- ☐ If it comes into contact with your eyes, rinse immediately with clear water and consult a doctor.

### Notes For Your Health

- ☐ Consult a dentist before using the product if you have any questions or doubts about the product and/or if you have sensitive teeth.
- ☐ Do not use the product if you have inflammation of the gums, gum lacerations or after undergoing a dental or jaw operation.
- ☐ If gum irritations or inflammations occur, consult a dentist immediately.

- ☐ Do not use the product if you are prone to allergic or sensitive reactions to the ingredients contained in the product.
- ☐ Teeth whiteners can cause increased sensitivity of the teeth (e.g. when consuming hot or cold drinks). Consult a dentist if you notice increased, persistent sensitivity.
- ☐ Do not use any additional cosmetic oral care products (e.g. whitening toothpastes or fluoride toothpastes etc.) in conjunction with the teeth whitening strips.

### Instructions for Use

- ☐ Do not allow the product to take effect for longer than 30 minutes.
- ☐ Do not apply more than once a day! There must be a gap of at least six hours between two applications.
- ☐ Each teeth whitening strip may only be used once.
- ☐ Do not use teeth whitening strips directly after brushing your teeth. There must be a gap of at least an hour between when you brush your teeth and when you apply the teeth whitening strips.
- ☐ Never expose the teeth whitening strips to direct sunlight or other heat sources.
- ☐ Store the teeth whitening strips in ambient temperatures of below 25°C.

## Items Supplied

- 14 x 2 teeth whitening strips

INCI: GLYCERIN; AQUA; CELLULOSE GUM; SODIUM CHLORITE; EDTA; D,L-MENTHOL; CITRIC ACID

## Use

- To achieve a better result, you are advised to restrict or avoid smoking and consuming coffee, tea, cola, red wine, cherries, blueberries and other highly pigmented foodstuffs.
  - Regular application every three months is recommended in order to maintain the effect permanently.
1. Peel off the curved teeth whitening strip from the plastic film.
  2. Place the teeth whitening strip with the gel side on the top teeth and carefully press it firmly on the teeth using your fingers. Make sure that the teeth whitening strip begins at the edge of the gum and covers the entire tooth.
  3. Repeat the first two steps with the straight teeth whitening strip for the bottom teeth.

4. Allow the teeth whitening strips to take effect for no more than 30 minutes.
5. Carefully peel the two teeth whitening strips off the teeth. If they are stuck on too firmly, moisten the teeth whitening strips with a little water.
6. Brush your teeth and rinse your mouth out thoroughly with clear water to remove any gel residues.

Note: If after using the teeth whitening strip white spots appear on your teeth, this is a temporary effect caused by uneven areas on the tooth. The white spots will disappear after a little while.

## Disposal

The packaging material is recyclable. Please dispose of it in an environmentally friendly manner and take it to a recycling centre. Dispose of the set in an environmentally friendly manner as well. You can obtain more information from your local authorities.

---

Customer Service / Importer:

DS Produkte GmbH, Am Heisterbusch 1, 19258 Gallin, Germany  
Tel.: +49 (0) 38851 314650 \*)

\*) Calls subject to a charge.  
All rights reserved.

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DS Produkte GmbH, Am Heisterbusch 1, 19258 Gallin, Duitsland